




NOTE:
ALL DIMENSIONS ARE APPROXIMATE
AND NOT INTENDED TO BE USED
FOR MANUFACTURE

STAND w/ 520mm TRACK
SCOPE HOLDER: Ø76mm
TRACK LENGTH: 520mm
TRAVEL DISTANCE: 385mm
STAND THROAT DEPTH: 120mm
COURSE FOCUS: 23.8mm/Rev
FINE FOCUS: 0.3 mm/Rev

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES [MM] AND INCLUDE CHEMICALLY APPLIED OR PLATED FINISHES		PROPRIETARY AND CONFIDENTIAL THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF NAVITAR. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF NAVITAR IS PROHIBITED.		 NAVITAR ROCHESTER, NEW YORK				
SCALE: 1:4				STAND, 76mm HOLDER, FINE FOCUS, HIGH TRACK				
WEIGHT: 8.27 lbmass				SIZE B	PART NUMBER 8-28256	REV		
DR BY: C. HATCH								
CH BY: DEROSA		DO NOT SCALE DRAWING						
PROJECT #: E3592								
DATE DR: 4/24/2018								

OUTLINE