



NA = 0.26  
WD = 30.5mm



## OUTLINE

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES AND INCLUDE CHEMICALLY APPLIED OR PLATED FINISHES		 ROCHESTER, NEW YORK	
SCALE:	1:1		
WEIGHT:	200 grams	OBJECTIVE, 10X, M-PLAN, APO, NIR	
DR:	R. PODLENA		
CH:		SIZE	REV
PROJECT #:		A	A
DATE DR:	5/16/2002	1-62238 SHEET 2 OF 2	

THE CONTENTS OF THIS  
DOCUMENT ARE PROPRIETARY  
TO NAVITAR, ITS USE IS  
AUTHORIZED ONLY FOR  
RESPONDING TO REQUEST FOR  
QUOTATION OR PERFORMANCE  
OF WORK FOR NAVITAR

0.0 inch 1.0  
0.0 mm 25.0

DO NOT SCALE DRAWING

